


## Technical data

### PART ONE

### Evaluation

#### 1. CT Value and TCID50 Conversion

<b>CT Value</b>	19.01	22.46	25.225	29.03	32.195	36.385	39.35
<b>TCID 50</b>	200,000	20,000	2,000	200	20	2	0.2



#### 2. SARS-CoV-2 virus evaluation procedure and results

2.1 **Location:** Biological safety protection third-level laboratory

2.2 **Virus:** SARS-CoV-2 WIV04

2.3 **Virus amplification cell:** Vero E6

2.4 **Titer detection:** limiting dilution assay, the origin virus' titer  $2 \times 10^6$ TCID<sub>50</sub>/ml

2.5 **Test Kit:** BALLYA Covid 19 Antigen Rapid Test Kit

2.6 **Virus dilution:**

2 x 10<sup>5</sup>TCID<sub>50</sub>/ml

2 x 10<sup>4</sup>TCID<sub>50</sub>/ml

2 x 10<sup>3</sup>TCID<sub>50</sub>/ml

2 x 10<sup>2</sup>TCID<sub>50</sub>/ml

2 x 10<sup>1</sup>TCID<sub>50</sub>/ml

2 TCID<sub>50</sub>/ml

2.7 **Sample add:** 80μl / cassette

2.8 **Detection time:** 15 minutes

2.9 **Repeat:** each sample repeat 3 times

**Results**

No.	Virus Concentration	Sample added volume (80µl / cassette)	Results
1	2 x 10 <sup>5</sup> TCID <sub>50</sub> /ml	1.6 x 10 <sup>4</sup> TCID <sub>50</sub> /ml	High Brightness line
2	2 x 10 <sup>4</sup> TCID <sub>50</sub> /ml	1.6 x 10 <sup>3</sup> TCID <sub>50</sub> /ml	High Brightness line
3	2,000 TCID <sub>50</sub> /ml	160 TCID <sub>50</sub>	Medium brightness line
4	200 TCID <sub>50</sub> /ml	16 TCID <sub>50</sub>	Weak line
5	100 TCID <sub>50</sub> /ml	8 TCID <sub>50</sub>	Weak line
6	20 TCID <sub>50</sub> /ml	1.6 TCID <sub>50</sub>	No line
7	2 TCID <sub>50</sub> /ml	0.16 TCID <sub>50</sub>	No line

**Conclusion:**

The sensitivity **8 TCID<sub>50</sub> (32.195 < CT value < 36.385)**

**PART TWO**

**National Reference & Company Standard**

**National reference (China)**

**Negative reference:** total 20

(N1,N2,N3,N4,N5,N6,N7,N8,N9,N10,N11,N12,N13,N14,N15,N16,N17,N18,N19,N20)

**Positive reference:** total 8 (P1,P2,P3,P4,P5,P6,P7,P8)

**Minimum detection reference:** total 1 (S)

**Repeat reference:** total 2 (R1, R2)

**Company Standard**

**BALLYA**

**1. Minimum detection : 7 TCID<sub>50</sub>/ml**

Use **National Reference Panel for 2019-nCoV Antigen Detection Kit (S)** ,

**1:140,000** dilution or titer 7 TCID<sub>50</sub>/ml as company reference, 20 times detection in same batch, positive results not less than 95%

**2. Positive rate:** 100% (detection P1-P8, all positive)

**3. Negative rate:** 100% (detection N1-N20)

**4. Repeat:** all positive results (R1= 10 times detection, R2= 10 times detection); color band of cassette in same brightness

**5. Batch CV%:** use three different batches to detect R1 & R2, each reference 10

- times detection, all positive results; color band of cassette in same brightness
6. **Stability:** Accelerated test, 55°C for 100 days (equivalent to); then detect according to No.1 - No.4; all results are match.
  7. **Hook effect:** dilute 1:1 national reference (S) or company reference 0.5x10<sup>6</sup> TCID<sub>50</sub>/ml, no hook effect

**Cross reactivity:**

There is no cross-reaction with pathogens at the following concentrations.

No.	Phathogen	Concentration
1	HCoV-HKU1	1.5x10 <sup>6</sup> copies/mL
2	HCoV-OC43	1.1x10 <sup>6</sup> copies/mL
3	HCoV-NL63	1.0x10 <sup>6</sup> copies/mL
4	HCoV-229E	3.8x10 <sup>6</sup> copies/mL
5	Novel influenza A(h1n1) virus(2009)	1.8x10 <sup>7</sup> copies/mL
6	Seasonal h1n1 influenza virus	8x10 <sup>7</sup> copies/mL
7	Influenza A virus(h3n2)	1.2x10 <sup>7</sup> copies/mL
8	Influenza A virus(h5n1)	10 <sup>6</sup> copies/mL
9	Influenza A virus(h7n9)	10 <sup>6</sup> copies/mL
10	Influenza B (yamagata)	2.8x10 <sup>7</sup> copies/mL
11	Influenza B (victoria)	2.0x10 <sup>7</sup> copies/mL
12	Respiratory syncytial virus type A	10 <sup>6</sup> copies/mL
13	Respiratory syncytial virus type B	1.2x10 <sup>6</sup> TCID <sub>50</sub> /mL
14	Parainfluenza virus type 1	10 <sup>6</sup> copies/mL
15	Parainfluenza virus type 2	10 <sup>6</sup> copies/mL
16	Parainfluenza virus type 3	10 <sup>6</sup> copies/mL
17	Rhinovirus A	10 <sup>6</sup> copies/mL
18	Rhinovirus B	10 <sup>6</sup> copies/mL
19	Rhinovirus C	10 <sup>6</sup> copies/mL
20	Adenovirus 1	10 <sup>6</sup> copies/mL
21	Adenovirus 2	10 <sup>6</sup> copies/mL
22	Adenovirus 3	10 <sup>6</sup> copies/mL
23	Adenovirus 4	10 <sup>6</sup> copies/mL
24	Adenovirus 5	10 <sup>6</sup> copies/mL
25	Adenovirus 7	10 <sup>6</sup> copies/mL
26	Adenovirus 55	10 <sup>6</sup> copies/mL
27	Human metapneumovirus	10 <sup>6</sup> copies/mL
28	Enterovirus A	10 <sup>6</sup> copies/mL
29	Enterovirus B	10 <sup>6</sup> copies/mL
30	Enterovirus C	10 <sup>6</sup> copies/mL
31	Enterovirus D	10 <sup>6</sup> copies/mL
32	EB virus	10 <sup>6</sup> copies/mL
33	Measls virus	10 <sup>6</sup> copies/mL

34	Human cytomegalovirus	10 <sup>6</sup> copies/mL
35	Rotavirus	10 <sup>6</sup> copies/mL
36	Norovirus	10 <sup>6</sup> copies/mL
37	Mumps virus	10 <sup>6</sup> copies/mL
38	Herpes zoster virus	10 <sup>6</sup> copies/mL
39	Mycoplasma pneumoniae	10 <sup>6</sup> CFU/mL

**Interfering substances:**

The kit is not affected with the substances at the following concentrations

No.	Interference substances	Concentration
1	Mucin	60mg/dL
2	Whole blood	20%(v/v)
3	Phenylephrine	2mg/mL
4	Oxymetazoline	2mg/mL
5	Sodium chloride (with preservative)	0.1%
6	Beclomethasone	20mg/mL
7	Dexamethasone	20mg/mL
8	Flunisolide	20μg/mL
9	Triamcinolone acetonide	2mg/mL
10	Budesonide 2mg/mL 23 Ceftriaxone	40μg/mL
11	Mometasone	2mg/mL
12	Fluticasone	2mg/mL
13	Zanamivir	20mg/mL
14	Peramivir	1mg/mL
15	Lopinavir	500mg/mL
16	Ritonavir	60mg/mL
17	Interferon-α	800IU/mL
18	Ribavirin	10mg/mL
19	Oseltamivir	60ng/mL
20	Arbidol	700ng/mL
21	Levofloxacin	10μg/mL
22	Azithromycin	1mg/L
23	Ceftriaxone	40μg/mL
24	Meropenem	200mg/mL
25	Tobramycin	0.6mg/mL
26	Histamine dihydrochloride	5mg/mL