





PIPETTE TIP SORTER

For fully automatic quick filling of pipette tips

6 Boxes Only 7 minutes



Pipette tip sorter is an exclusive innovative and patented product, designed to meet the needs of biological laboratories. It can automatically dispense pipette tips into tip boxes, which is an efficient solution.



3 points:

Fully automated

Suitable for 2 pipette tip types

more environmentally friendly

7 advantage:

- ① No vibration at all, no influence on other instruments;
- ② 10 million boxes of continuous pipette tip experiment without calibration;
- ③ low noise;
- ④ 2 kinds of tips can be installed without replacing parts;
- ⑤ Can be placed on the desktop, 67*56*51cm, 1.5 microwave oven size;
- ⑥ Can prevent jam, 19+ sensors, anti-stacking, anti-double insertion and anti-extrusion;
- ⑦ Modular quick-release feeder tube, users can easily maintain.

parameter:

Size	W 670mm D 560mm H 510mm
Weight	40KG
Material	Aviation grade aluminum alloy + medical grade 304 stainless steel
power input	220V-12V 7A
Operating power	40W \pm 15%
Support tip	Axygen, Eppendorf and other mainstream brands
Platform move	X:560mm y:540mm
Storage bin volume	200ul:1500 10ul:2500
Ambient temperature	15-30°C
Relative humidity	<80%RH
Filling time	0.8 seconds/1 tip

15 characteristics:



Faster speed



Stable and long-lasting operation



Aerospace-grade aluminum alloy



German quality



Graphical operation interface



Trinamic drive



Industrial-grade power supply



Full brushless system



Non-contact sensing



Sterilizable inner warehouse



Mute function of the whole machine



Magnetic suction platform



German quality



Quickly dismantle the material tube



Larger capacity

BALLYA

Address: 5F, Building B, No.241, ShangYe Rd., Huadu, Guangzhou, China

Sales: +86 20 3947 9163

Email: info@ballyabio.com

Web: <https://ballyabio.com/>